# ADA S7/IDS

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# ADA S7/IDS

IDS pH adapter for pH electrodes  
with S7 plug-in connector

**Operating manual**

**Note**

The latest version of the present operating manual can be found on the Internet under [www.WTW.com](http://www.WTW.com).

**Copyright**

© Weilheim 2010, WTW GmbH

Reprinting - even as excerpts - is only allowed with the explicit written authorization of WTW GmbH, Weilheim.

Printed in Germany.

## General information

The ADA S7/IDS adapter forms an IDS pH sensor in conjunction with a pH electrode with S7 plug-in connector.

### Automatic adapter recognition

The measuring electronics with the stored adapter data is in the adapter head of the ADA S7/IDS. The data include, among other things, the adapter type and series number. With each calibration, the calibration data is written in the adapter and the calibration history is recorded. The data is recalled by the meter when the adapter is connected and is used for measurement and for measured value documentation. Storing the calibration data in the adapter ensures that the correct slope and asymmetry are automatically used if the sensor is operated with several meters.

The digital transmission technique guarantees the failure-free communication with the meter even with long connection cables. If the sensor firmware is enhanced by WTW, it can be updated via the meter.

#### Note

More detailed information on the automatic adapter recognition or sensor recognition is given in the operating manual of the meter.

### Calibration data and calibration history

The connected meter uses for measurement the calibration data stored in the adapter, regardless of the mounted electrode.

#### Note

Therefore, you have to calibrate the IDS pH sensor each time the electrode with the S7 plug-in connector has been exchanged. Otherwise, the IDS pH sensor will work with the calibration data of the previously used electrode. This could lead to incorrect measured values. When evaluating the calibration history, please note that in between calibrations, the electrode might have been exchanged.


## Technical data

<b>General data</b>	Electrodes that can be connected		All electrodes with S7 plug-in connector	
<b>Application characteristics</b>	pH measuring range		Depending on the used electrode	
	Allowed temperature range		0 ... 80 °C (32 ... 176 °F)	
	Typical application		Applications with special pH electrodes	
<b>Characteristics of the IDS measuring technique</b>	<b>Measured parameter</b>	<b>Measuring range</b>	<b>Resolution</b>	<b>Accuracy (± 1 digit)</b>
	pH	-2.000 ... 20.000	0.001	± 0.004
	U [mV]	-1000.0 ... +1000.0	0.1	± 0.2
<b>Material, electrical connection</b>	Material of the plug enclosure		PPE/PS	
	Combination electrode connection		S7 plug	
	Meter connection		Fixed cable with digital plug	
<b>Connection cable</b>	Length		1.5 m	
	Diameter		4.3 mm	
	Smallest allowed bend radius		Fixed installation:	20 mm
			Flexible use:	60 mm
<b>Electrical data of the electrode connection</b>	Plug type		Socket, 4 pins	
	Input resistance		> 5 * 10 <sup>12</sup> ohm	
	Input current		< 1 * 10 <sup>-12</sup> A	

## Commissioning, measuring, calibration

### Commissioning

Connect the S7 plug-in connector of the pH electrode with the S7 plug of the adapter.

#### Note

The further commissioning steps are given in the operating manual of the pH electrode.

- Connect the electrode to the meter.
- Calibrate the electrode according to the operating manual of the meter and observe the following rules while doing so:

### Calibration and measurement

Connect the adapter with the pH electrode to the meter. All further details are given in the operating manuals of the pH electrode and meter.

**Caution: Heed the depth of immersion of the pH electrode. The ADA S7/IDS adapter must not be submersed!**

## Storage

For storage, leave the adapter plugged on the pH electrode if you do not intend to use both devices otherwise. The storing method depends on the pH electrode. More details are given in the operating manual of the pH electrode.

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## **Wissenschaftlich-Technische Werkstätten GmbH**