

Keofitt A/S

Hillerød, 17th september 2018

Kullinggade 31 B+E

5700 Svendborg

Certificate

We hereby confirm that all Keofitt products mentioned have undergone the surface treatment described in the table below.

Surface treatment	1. Electropolishing 2. Decontamination 3. Rinsing under pressure
Chemicals applied	1. Poligrat E269 for electropolishing 2. HNO₃ for decontamination 3. Distilled water (pure water) for rinsing
Mounting	Items are mounted in fixtures made from titanium
Treatment time	1. 2 to 7 minutes for electro polishing depending on material quality 2. Minimum 30 minutes for decontamination 3. Minimum 2 minutes for rinsing
Article number	Apply to all Keofitt valves and accessories made from steel and according to drawings and specifications.

Morten Clausen, Managing director