





PP BOTTLE 2000 ml


ART. NO. 260625

GENERAL

-  KEOFITT has the widest selection of spare parts and accessories to complete your sampling system.
-  Bottles and closures are included in parts & accessories - contact KEOFIT for further details.
-  Autoclavable.
Recommended autoclave cycle is 121°C, 15 psig (1 bar) for 20 minutes. For best results, use a slow exhaust cycle.
-  When autoclaving, the bottle cannot be sealed. Remember to loosen or remove the closure.



FEATURES

-  Closure included (#260328)

CERTIFICATION*

· FDA CFR 21 §177.1520 · USP Class VI

TECHNICAL DATA

Material (process contact)

Material (without process contact)


- Closure: Polypropylene (PP)
- Gasket: Thermoplastic elastomer (TPE)
- Bottle: Polypropylene copolymer (PPCO)
- Temperature range: -20° - +121°C / -4° - +249.8° F
- Pressure: Sampling should be at zero pressure
- Microwaveability: Marginal
- Sterilization:
 - Autoclaving: Yes
 - EtO gas: Yes
 - Dry heat: No
 - Radiation: No
 - Disinfectants: Yes
- Chemical Resistance:
 - Not recommended for: Aromatic ketones - Ethers - Aromatic/Halogenated Hydrocarbons

Net Weight

- kg/lbs: 0.313 kg/0.69 lbs

Spareparts

- 260328: CLOSURE f/PP BOTTLE 250/1000/2000 ml
- 260327: 10 PACK GASKET f/CLOSURE f/PP BOTTLE 250/1000/2000 ml



*For further information and download please visit www.keofitt.dk

KEOFITT A/S reserves the right to change technical data without notice! Valid version of the data sheets is available on www.keofitt.dk.