





# PP BOTTLE 4000 ml


ART. NO. 260725

## GENERAL

-  KEOFITT has the widest selection of spare parts and accessories to complete your sampling system.
-  Bottles and closures are included in parts & accessories - contact KEOFITT for further details.
-  Autoclavable.  
Recommended autoclave cycle is 121°C, 15 psig (1 bar) for 20 minutes. For best results, use a slow exhaust cycle.
-  When autoclaving, the bottle cannot be sealed. Remember to loosen or remove the closure.



## FEATURES

-  Closure included (#260428)

## CERTIFICATION\*

· FDA CFR 21 §177.1520 · USP Class VI

## TECHNICAL DATA

### Material (process contact)

### Material (without process contact)


- Closure: Polypropylene (PP)
- Gasket: Thermoplastic elastomer (TPE)
- Bottle: Polypropylene copolymer (PPCO)
- Temperature range: -20° - +121°C / -4° - +249.8° F
- Pressure: Sampling should be at zero pressure
- Microwaveability: Marginal
- Sterilization:
  - Autoclaving: Yes
  - EtO gas: Yes
  - Dry heat: No
  - Radiation: No
  - Disinfectants: Yes
- Chemical Resistance:
  - Not recommended for: Aromatic ketones - Ethers - Aromatic/Halogenated Hydrocarbons

### Net Weight

- kg/lbs: 0.621 kg/1.368 lbs

### Spareparts

- 260428: CLOSURE f/PP BOTTLE 4000/5000 ml
- 260427: 10 PACK GASKET f/CLOSURE f/PP BOTTLE 4000/5000 ml



\*For further information and download please visit [www.keofitt.dk](http://www.keofitt.dk)