



# **VALVE BODY 4KAI TYPE MINI TRI-CLAMP**

ART. NO. 440021

### **GENERAL**



The KEOFITT 4KAI Sampling Valve is our low cost sterilizable sampling valve for less demanding applications. Replaces non-hygienic sampling cocks. More than 20 different configurations available.



The sampling valve can be used for any process sampling for chemical and/or physical analysis.



Cleaning/sterilizing: Between batches: Valve in open position: Cleanable by means of CIP using the detergent solution suitable for the actual process media. Between samples: Valve in its normal closed position: cleanable by CIP as "Between batches" or the valve may be sterilized by means of steam SIP or chemical SIP using a procedure appropriate to the actual circumstances. For further advice, please contact Keofitt.



Designed for sampling of liquids with a viscosity of up to approx. 100 cP containing no particles larger than  $\emptyset$ 1.5 mm. Sampling of more viscous liquids is possible, only will it take longer (depending on process pressure).



### **FEATURES**



Installation: Mini Tri-Clamp DIN 32676:2009-05 Table 4 Row C



Operation: Depending on choice of valve head



Membrane: Depending on choice of valve head

59 mm

Ø25 mm



Inlet / Outlet: Hose piece

## **CERTIFICATION\***

- · EU EC 1935/2004 · EU EC 2023/2006 · DK No. 1248 10/2018 · 3-A Certificate · 3.1 Material Certificate
- · ATEX 2014/34/EU · PED 2014/68/EU · Keofitt DoC

### **TECHNICAL DATA**

## **Material (product contact)**

· Steel parts AISI 316L (1.4404)

# **Material (without product contact)**

Steel parts
Rubber cap
Stainless steel chain
AISI 316L (1.4404)
EPDM (black)
AISI 316L (1.4404)

### **Surface Treatment**

 $\begin{array}{ll} \cdot \text{Outside} & \text{Ra} \leftrightharpoons 1.2 \, \mu\text{m} \\ \cdot \text{Inside (wetted surface)} & \text{Ra} \leftrightharpoons 0.8 \, \mu\text{m} \\ \cdot \text{Process connection} & \text{Ra} \leftrightharpoons 0.8 \, \mu\text{m} \end{array}$ 

#### **Pressure & Temperature**


Pressure
Temperature
Air supply
Depending on choice of valve head
Depending on choice of valve head

#### **Net Weight**

· kg/lbs 0.15 kg/0.331 lbs

## **Spareparts**

774401 PARTS FOR 4KAI VALVE BODY







38