



VALVE HEAD W9 TYPE H

ART. NO. 600041

GENERAL



The KEOFITT CLASSIC W9 Sampling Valve is the original and leading sterilizable sampling valve in the world. Used in all industries for decades. More than 320 standard valve configurations. Unique serial no. for each valve*.



The sampling valve can be used for any process sampling for microbiological, chemical and/or physical analysis.



Cleaning/sterilizing: Between batches: Valve in open position: Cleanable by means of CIP using the detergent solution suitable for the actual process media. Between samples: Valve in its normal closed position: cleanable by CIP as "Between batches" or the valve may be sterilized by means of steam SIP or chemical SIP using a procedure appropriate to the actual circumstances. For further advice, please contact Keofitt.

Not recommended for autoclave due to plastic parts.



Designed for sampling of liquids with a viscosity of up to approx. 1.000 cP containing no particles larger than Ø3 mm. Sampling of more viscous liquids is possible, only will it take longer (depending on process pressure).



FEATURES



Installation: Threaded socket connection M22x1.5



Operation: Turn knob - spring loaded (opens clockwise)



Membrane: Silicone (#600051)

CERTIFICATION*

- · EU EC 1935/2004 · EU EC 2023/2006 · DK No. 1248 10/2018 · 3-A Certificate · ATEX 2014/34/EU
- · PED 2014/68/EU · FDA CFR 21 §177.2600 · USP Class VI · Keofitt DoC

TECHNICAL DATA

Material (process contact)

· Membrane SILICONE, GREY (#600051)

Material (without process contact)

· Steel parts AISI 303 (1.4305) / AISI 316L (1.4404)

· Knob PA6 (Black)

Membrane

· Article SILICONE, GREY (#600051)

· Certification* · EU EC 1935/2004 · EU EC 2023/2006

 \cdot DK No. 1248 10/2018 \cdot FDA CFR 21 §177.2600 \cdot USP Class VI

· REACH · RoHS · ADI Free · Keofitt DoC

Pressure & Temperature

• Pressure 0 - 10 bar / 0 - 145 psi • Temperature 1 - 130°C / 34 - 266°F • Air supply


Net Weight

· kg/lbs 0.228 kg/0.501 lbs

Spareparts

776041 PARTS FOR W9 HEAD 600041/42/43/47







^{*}For further information and download please visit www.keofitt.dk

KEOFITT CLASSIC W9 SAMPLING VALVE - BODIES

| KEC | | | M9 BODIES | | | | ON RE |
|--|-------------------|---------------------|--|------------------|-----------------|-------------------------------------|------------------|
| OFII | | | HOSE PIECE | 7\2" TRI-CLAMP | сүтхөтм дазянт | T\S. METDING ENDS | ON REQUEST ONLY! |
| TT CLASSIC | | TANK WELDING Ø28 | 100058 | 850008 | O STOOSE STOOSE | 850007 | NLY! |
| KEOFITT CLASSIC W9 SAMPLING VALVE - BODIES | | PIPE WELDING 1" | 120028 | 85010 | O Spoots | 850019 | |
| | WELDING BODIES | PIPE WELDING NW25 | A20028 | 15 anto Juo | 9100SB | ¹⁵ anba _l (10 | -34- |
| | | IN-LINE VERTICAL | 880027 | 15 antos 40 | 15 an to 140 | ISanbay UO | |
| | | IN-LINE HORIZONTAL | esoose. | 15 anto 140 | 15 an to 140 | ISanbas UO | <u>v</u> |
| | | | | | | | |
| | | VARIVENT Ø50 | O RESCOOSE | 85020 | Isanta, uo | Isanbar uo | \(\bar{A} \) |
| | | VARIVENT @68 | 6000SB | 85013 | 15311to140 | ISANDAJ UO | |
| | CONNECTION BODIES | MINI CLAMP | SECONS. | 35039 | Isanbaruo | Isanba, uo | |
| | N BODIES | CLAMP 1" & 116" | C 120059 | 850022 | 050024 | 15anbar uo | ŚŹ |
| | | CLAMP 2" | 0 S S S S S S S S S S S S S S S S S S S | 850026 850026 | 15anbal UO | Isanbar uo | |
| | | CLAMP 3" | D0008 | 15anto 140 | 15anta140 | sanba, uo | |
| | | TREAD | M28/1.5 850031 BIN 11651 DN 10 650040 | | M28XLS | | |
| | | | | | AND AND | | |