

VALVE BODY W9 TYPE P W/TC



ART. NO. 850010

GENERAL



The KEOFITT CLASSIC W9 Sampling Valve is the original and leading sterilizable sampling valve in the world. Used in all industries for decades. More than 320 standard valve configurations. Unique serial no. for each valve*.

The sampling valve can be used for any process sampling for microbiological, chemical and/or physical analysis.



Cleaning/sterilizing: Between batches: Valve in open position: Cleanable by means of CIP using the detergent solution suitable for the actual process media. Between samples: Valve in its normal closed position: cleanable by CIP as "Between batches" or the valve may be sterilized by means of steam SIP or chemical SIP using a procedure appropriate to the actual circumstances. For further advice, please contact Keofitt.



Designed for sampling of liquids with a viscosity of up to approx. 1.000 cP containing no particles larger than Ø3 mm. Sampling of more viscous liquids is possible, only will it take longer (depending on process pressure).

FEATURES



Installation: Pipe welding: DIN EN ISO 1127 25x1.25



Inlet / Outlet: 1/2" Clamp DIN 32676:2009-05 Table 4 Row C



· EU EC 1935/2004 · EU EC 2023/2006 · DK No. 1248 10/2018 · 3-A Certificate · 3.1 Material Certificate · Ra Certificate · ATEX 2014/34/EU · PED 2014/68/EU · Keofitt DoC

TECHNICAL DATA

Material (process contact)

· Steel parts

Material (without process contact)

 \cdot Steel parts

AISI 316L (1.4435 BNII) / AISI 316L (1.4404)

Surface Treatment

· Outside

Inside (wetted surface)
Process connection

Electropolished Ra <= $1.2 \,\mu$ m Electropolished Ra <= $0.5 \,\mu$ m Electropolished Ra <= $0.5 \,\mu$ m

Depending on choice of valve head

Depending on choice of valve head

Depending on choice of valve head

AISI 316L (1.4435 BNII)


Pressure & Temperature

- Pressure
- Temperature
- Air supply
- **Net Weight**
 - ∙ kg/lbs

Spareparts


0.388 kg/0.855 lbs





Operation: Depending on choice of valve head

Membrane: Depending on choice of valve head



*For further information and download please visit www.keofitt.dk

KEOFITT A/S reserves the right to change technical data without notice! Valid version of the data sheets is available on www.keofitt.dk.

Last updated 01-05-2020

KEOFITT CLASSIC W9 SAMPLING VALVE - HEADS

Type B No spring 0-12 bar(g)				
Type N Preumatic, Adjustable 0-6 bar				
Type N Procumatic 0-6 bar				
Type Q Lever handle	C00003			
Type K Kay ring with Microport option	600049		DIQ OT OND JON	old Blord Joh
Type K Key ring			470000 AZ PO000	§
Turn knob with Microport option			olden ever john	OIGBIELD JOH
Type H Turn knob	1 70000		Throom	
	sincolli2	EPDM	EEKW	PTFE