


PRODUCT CATALOG 2021

KEOFITT MICRO PORT



KEOFITT®
WORLD LEADERS IN STERILE SAMPLING™

SAMPLING



KEOFITT MICRO PORT

The KEOFITT MICRO PORT is a re-sealing rubber septum designed to take samples using hypodermic needles or syringes with minimal risk of contamination. Designed for sampling low viscosity products in the range 0-50 cP.

The KEOFITT MICRO PORT is easy to use and maintain, if needle sampling is required.






The Micro Ports are used in a range of business areas, such as brewery, dairy, pharmaceutical and biotechnological industries, with a wide range of products such as:




- Water
- Medical solutions
- Milk
- Filtered Beer
- Perfume
- Olive Oil

Available as a single septum or as a MULTI MICRO PORT with seven septums.























KEOFITT MICRO PORT

MICRO PORT			
	TANK	PIPE WELDING 1”	PIPE WELDING NW25
MICRO PORT / MULTI MICRO PORT 49	 900056	 900057	 900014
	MINI CLAMP	CLAMP 1” & 1½”	VARIVENT
	 900059	 900058	On request

MULTI MICRO PORT 49		
TANK WELDING 7XBUTYL MEMBRANE	VARIVENT Ø68 7XBUTYL MEMBRANE	CLAMP 3” 7XBUTYL MEMBRANE
 840001	 840009	 840020

KEOFITT MICRO PORT - PARTS & ACCESSORIES

MICRO PORT PARTS & ACCESSORIES	O-RINGS	 900822 EPDM 9.3X2.4	 900824 EPDM 14X4	 900825 EPDM 44.6X2.4	 400208 EPDM 60X3				
	SEPTUM	 900049 BUTYL							
	MISC.	 840055 FILM MULTI	 840118 TOOL MULTI	 900054 NEEDLE	 900055 WICK	 900077 CHAIN	 600065 CHAIN FOR MULTI MICRO PORT	 900075 CLAMP RING 1/2"	 900086 CLAMP RING 1"
		 900074 GASKET EPDM 3/4	 900087 GASKET EPDM 1/2"	 900091 GASKET EPDM 1"	 900093 GASKET 3" CLAMP EPDM	 840102 INNER PLATE FOR MICRO PORT	 900102 INSERT FOR MICRO PORT	 900441 PLUG	